

Shimin config file

DataPath: #Not sure if open to users

proj_path_base : "/SNS/SNAP/IPTS-26894/shared/autoreduce/mcp" # Path where save project data in MCP

button to browse for input folder

OB_path_base : "/SNS/SNAP/IPTS-26894/shared/autoreduce/mcp/OB_May2022" # Path where save OB data in MCP

button to browse for OB folder

result_path : "/SNS/SNAP/IPTS-26894/shared/insitu_recon/output_loopCronTest_10C_Tilt8" # Path where save the output data

~~output folder or~~ provide the location for them.

search_string: "Tilt8deg_CT_GS_loop_10C" #SHOULD BE SAME AS API folder_tag

same as folder_tag, should we figure it out, provide a default

log_name : '/SNS/SNAP/IPTS-26894/shared/insitu_recon/logs/reclog_loopcrontest.log' # Not sure

~~—# provide option in advanced tab to customize the location and name of the log file name~~

""OB Parameter (May open to advanced users **NOT SURE)""

OB_sub_folder_exists : False #open beam data path sub-folders

what is that - go and look inside the sub-folders.

num_OB_set : 5 #number of open beam datasets

number of ob to use

in advanced tab for sure. It will be a radio-button

OB_gain_factor : 1.00 #gain factor for open beam to match the counts of projection data

in advanced. Should we provide an option to calculate this one by doing ROI on integrated signal

auto-populate this by finding the proton charge of sample and ob

""""""

#SNAP API params

SnapApi: # Need be set at the beginning of measurement to start experiment manually and saved for following code running

goni1: 8.0 # degrees

what is the units. Which one is goni1. Should we provide a picture? This should go in a second tab called "instrument"

goni2: 0.0 # degrees

same thing

folder_tag: 'Tilt8deg_CT_GS_loop_10C' # "IMPORTANT" same as searching string

no need to have it twice in the UI

file_tag: '8degtilt'. # not important but keep in advanced tab (temporary file used in reconstruction)

how is the tag used. Can we figure it ourselves?

run_title: 'Tilt CT' # display in the queue on SNAP interface , name of the job

```

# where and how this is used
run_notes: '8 degree tilts at 10 C'. # same thing as above
# probably with the first tab. TBD
p_charge: 10.0 # maybe we can try to retrieve it from the sample data
# in the tab "instrument"

# Evaluation ROI

roi_info : [80, 100, 200, 200] # Should be determined by users but need visualization setting
after 1st reconstruction
[ymin, xmin, height, width]
# the UI should have a button saying RUN FIRST RECONSTRUCTION as this option
should be enabled after the first reconstruction is done
# if None, take everything

Qlthresh: 4.3 # User Setting
# what is that

# in experiment tab
WaveOpts: # User setting
wav_idx_start : [1230]
wav_idx_end : [1271]

Detector offset
Distance source detector
# Should we instead ask for the range in lambda or TOF units and figure out ourselves that
index range
# that means we will need to add the Detector offset and distance source/detector

# in advanced tab
BasicPara:
detector_rows : 512 #May from Meta data, open to advanced users
detector_columns : 512 #May from Meta data, open to advanced users

# cropping
z_start : 150 #User setting # what is that. How it's used. # start
z_numSlice : 250 #User setting # how it's used. # how many to go
cut_cols : 400 #User setting # how it's used, how to label it. # width to take from center
(maybe replaced by next row)
# maybe new needed. Col_start and col_num

# be able to display profile (counts vs lambda/tof/file index)
# be able to display projections

# should we hard code those?
# if not, what is their meanings?

```

advanced for imaging team only

align_paras: # Hide to users

x_off : 2

y_off : 2

cent_x : 256

cent_y: 256

unused_pix_num : 1

neighbor_pix_num : 3

advanced?

more detail about what they are in case we need to change them

mbir parameters

reconstruction parameters

RecOpts: # Hide to users

rot_center: 254.25 # get this one from the UI using 0 and 180degrees (just like pyMBIR_UI)

T : 2.0

p : 1.2

sharpness : 0.0

snr_db : 30.0

max_iterations : 200

num_threads : 120

not mbir parameters (move to section that can be changed every time)

what is that?

Max_ang_num: 55 . # max number of measurement

what is that

FreqStop: 3. # number of acquisition before running evaluation code again